

αSyn does not impact *Enterococcus faecalis* planktonic growth




Fig 1. *E. faecalis* planktonic growth curve: (A) *E. faecalis* planktonic growth measured over time in the presence of varying [αSyn]. Conditions: Brain heart infusion media (BHI), 37 °C, (n = 5 biological, 3 technical replicates; error bars are + SEM). Data fitted to Gompertz growth eqn using non-linear regression, straight line (B), (C), and (D). Parameters for Gompertz determined from fits. Data are shown as mean ± std dev. Data compared using 2-way ANOVA (** = p < 0.005).

- *E. faecalis* growth was monitored over time in the presence of varying αSyn concentrations
- αSyn did not have a biologically significant impact on lag time, growth rate, doubling time, or maximum bacteria population

αSyn does not impact *Streptococcus mutans* planktonic growth




Fig 2. *S. mutans* planktonic growth curve: (A) *S. mutans* planktonic growth measured over time in the presence of varying [αSyn]. Conditions: BHI, 37 °C, (n = 5 biological, 1 technical replicate; error bars are + SEM). Data fitted to logistic growth eqn using non-linear regression, straight line (B), (C), and (D). Parameters for logistic determined from fits. Data are shown as mean ± std dev. Data compared using Welch's t test (** = p < 0.005).

- *S. mutans* growth was monitored over time in the presence of varying αSyn concentrations

• ~1 million people live with Parkinson's disease (PD) in the US. There are no disease-modifying treatments available.

• αSyn is associated with PD disease pathologies.

• Since evidence suggests gut and oral microbiome health influences PD pathology and αSyn is an antimicrobial agent, we examined the interaction between αSyn and gut/oral bacteria.

• Using methods to measure both planktonic and biofilm growth, we found αSyn differentially impacts gram positive bacteria.

• Contrary to previous results, we show αSyn can increase membrane integrity and biofilm formation in bacteria.

αSyn differentially affects *E. faecalis* and *S. mutans* biofilm formation




Fig 3. Biofilm Formation: (A) *E. faecalis* biofilm formation in BHI media with varying concentrations of α-synuclein following incubation at 37 °C for 48 hours. N = 5 biological replicates; 2 technical replicates. Error bars are ± SEM. (** significant at the level of p < 0.05) (B) *S. mutans* biofilm formation in BHI media with varying concentrations of α-synuclein following incubation at 37 °C for 48 hours. N = 5 biological replicates; 1 technical replicate. Error bars are ± SEM.

- Biofilm formation was evaluated following incubation for 48 hours in varying αSyn concentrations.
- αSyn increased the biofilm formation in *E. faecalis*, but had no significant impact on biofilm formation in *S. mutans*.

αSyn alters membrane integrity of *E. faecalis*




Fig 4. *E. faecalis* SYTO9/Propidium Iodide Fluorescence Assay: (A) Normalized *E. faecalis* SYTO 9/PI fluorescence at varied concentrations of α-synuclein. Data fit to one site binding – total and nonspecific binding (straight line) with 95% confidence interval (dashed line). N = 1 biological replicate with 3 technical replicates. (B) *E. faecalis* emission spectrum after addition of SYTO 9 and Propidium Iodide (PI) into varying ratios of untreated cells:cells treated with isopropanol. (C) Standard curve of SYTO 9/PI with percent of untreated cells

- SYTO 9 and Propidium Iodide (PI) fluorescence was observed at varying αSyn concentrations and normalized to fluorescence of increasing ratio of isopropanol treated cells.
- Increased SYTO 9/PI fluorescence ratio indicates increased membrane integrity of *E. faecalis* when exposed to α-synuclein

References

1. P. Baker, P. J. Hill, B. D. Snarr, N. Alnabehsaya, M. J. Pestrak, M. J. Lee, L. K. Jennings, J. Tam, R. A. Melnyk, M. R. Parsek, D. C. Sheppard, D. J. Wozniak, P. L. Howell, Exopolysaccharide biosynthetic glycoside hydrolases can be utilized to disrupt and prevent *Pseudomonas aeruginosa* biofilms. *Sci. Adv.* 2, e1501632 (2016).
2. G. A. O'Toole, Microtiter Dish Biofilm Formation Assay. *J. Vis. Exp.* (2011), doi:10.3791/2437.
3. Live/dead bacLight bacterial viability kits - thermo fisher scientific (2004), (available at <https://tools.thermofisher.com/content/sfs/manuals/mp07007.pdf>)

FUTURE DIRECTIONS

- We observed that α-synuclein is not bactericidal to the Gram-positive bacteria *Enterococcus faecalis* and *Streptococcus mutans* and our SYTO9/PI fluorescence assay demonstrated a positive correlation between α-synuclein concentration and increased SYTO 9:PI fluorescence, indicating enhanced bacterial viability.
- Given the relevance of α-synuclein to Parkinson's disease pathogenesis, future studies could explore how α-synuclein's effects on bacterial membrane integrity and biofilm formation alter host-microbe interactions and contribute to neurodegenerative processes in the context of Parkinson's disease

Acknowledgements

Gratitude to Tonya N. Zeczycki, Ph.D for providing me this opportunity to work in her lab. The experience I have gained is invaluable to my future medical career. Thanks to Whitney G. Bond for hands-on teaching and Mithil Patel for guidance. We thank the Wooten Brain Health Initiative Grant to TNZ and the Summer Scholar Research Program.